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**REF. 472 024**

## **I - GENERALITES**

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**Ordering Code: 472 024**

## **I – GENERAL**

This sensor is used to measure the instantaneous value of currents in the range of  $\pm 200 \mu\text{A}$  to  $\pm 100 \text{ mA}$  or the root mean square value of the sine component of currents in the range of 0 to 100 mA using the VTT system.

The measurements are performed via an input with safety banana jacks.

## **II – OPERATION**

### **1. VTT system**

This sensor is used in conjunction with the VTT console.

### **2. Sensor connection onto the VTT console**

The sensor is connected by means of the connector at the rear of the housing:

- Align the sensor with the selected measurement channel.
- Slide the sensor onto the guiding system in the console.
- Follow the guiding system of the VTT console.
- Push fully home.

### **3. Range**

The following ranges are available depending on the selector position:

- Measuring range from  $-100$  to  $+100 \text{ mA}$  (instantaneous value).
- Measuring range from 0 to  $+100 \text{ mA}$  (root mean square value).

### **4. Measurement type**

- Instantaneous value of the current.
- rms value of the sine component of the current.


## 5. Preparation for operation

- Connect the sensor.
- Select the type of measurement (instantaneous or rms) by means of the selector.
- Power up VTT.

## 6. Zeroing

If the zero is found to drift in either or both measurement types (instantaneous or rms), proceed as follows:

- Connect the sensor onto the console.
- Select the measurement type for which a drift is found.
- Power up VTT.
- On the sensor front, identify the access hole of the trim pot corresponding to the selected measurement type (instantaneous value / rms value).
- Use the screwdriver provided with the console to adjust the displayed value as follows:
  - $0 \pm 0.2$  mA for instantaneous current measurements.
  - $0 (+ 0.1 / -0.00)$  mA for rms current measurements.

## III – TECHNICAL DATA

- |   |                         |
|---|-------------------------|
| - Measuring range, instantaneous current: | from -100 mA to +100 mA |
| - Measuring range, rms current:           | from 0 mA to +100 mA    |
| - Resolution, instantaneous measurement:  | 0.2 mA                  |
| - Resolution, rms measurement:            | 0.1 mA                  |
| - Accuracy, instantaneous measurement:    | 1 % $\pm$ 0.4 mA        |
| - Accuracy, rms measurement:              | 1 % $\pm$ 0.2 mA        |
| - Bandwidth, instantaneous measurement:   | from 0 to 8 kHz         |
| - Bandwidth, rms measurement:             | from 30 Hz to 8 kHz     |
| - Input impedance:                        | less than 10 $\Omega$   |
| - Protection:                             | up to 1 A               |
| - Identification code:                    | 17                      |

## IV – MAINTENANCE

For any repair, adjustment or spare parts, please contact:

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S.A.V. JEULIN