Méthode VTT[®]

VTT[®] method

**Ref :
472 035**

Français – p 1

English – p 3

Version : 9007

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1 General

This sensor is used to measure the conductivity of solutions from $2 \mu\text{S}\cdot\text{cm}^{-1}$ up to $20 \text{mS}\cdot\text{cm}^{-1}$ by using the VTT system.

2 Operation

2.1 VTT system

This sensor is used in conjunction with the VTT console and with CTA probe.

2.2 Sensor connection onto the VTT console

The sensor is connected by means of the connector at the rear of the housing:

- Align the sensor with the selected measurement channel.
- Slide the sensor onto the guiding system in the console.
- Follow the guiding system of the VTT console.
- Push fully home.

2.3 Scale

Depending on the range selector setting, the following scales are available:

- from 0 to $2 \text{mS}\cdot\text{cm}^{-1}$ ($2 \text{mS}\cdot\text{cm}^{-1}$ range).
- from 0 to $20 \text{mS}\cdot\text{cm}^{-1}$ ($20 \text{mS}\cdot\text{cm}^{-1}$ range).

2.4 Measurement type

- Measurement of the electrical conductivity of aqueous solutions.

2.5 Preparation for operation

- Connect the CTA probe (ref. 701330) to the sensor.
- Connect the sensor onto the VTT console.
- Power up the VTT console.
- Select the desired range.
- Dip the probe into the solution so that the air escape holes are covered.
- Remove any air bubbles by softly stirring the probe in the solution.
- Wait for thermal stabilization of the probe (stable value on display).

2.6 Calibration

For correct measurement accuracy, the sensor must be calibrated prior to each experimental session.

2.6.1 2 mS.cm⁻¹ range

- Calibration is performed by using a 1413 $\mu\text{S.cm}^{-1}$ standard solution.
- Perform a measurement in the 2 mS.cm⁻¹ range by following the steps described in paragraph 5.
- Set the calibration button () until the display reads (1.413 \pm 0.004) mS.cm⁻¹.

2.6.2 20 mS.cm⁻¹ range

- Proceed as before using a standard solution to 12880 $\mu\text{S cm}^{-1}$ in order to succeed an accurate measure.

3 Technical data

- Measuring range (2 mS.cm⁻¹ range) : de 0 to 2 mS.cm⁻¹
- Measuring range (20 mS.cm⁻¹ range) : de 0 to 20 mS.cm⁻¹
- Resolution (2 mS.cm⁻¹ range) : 2 $\mu\text{S.cm}^{-1}$
- Resolution (20 mS range) : 20 $\mu\text{S.cm}^{-1}$
- Accuracy at 25°C (2 mS.cm⁻¹ range) : 3 % \pm 4 $\mu\text{S.cm}^{-1}$
- Accuracy at 25°C (20 mS.cm⁻¹ range) : 3 % \pm 40 $\mu\text{S.cm}^{-1}$
- Thermal drift between 15 and 35°C : \pm 0.5% /°C
- Response time from 20 to 35°C) : approx. 4 minutes
- Identification code (2 mS.cm⁻¹ range) : 33
- Identification code (20 mS.cm⁻¹ range) : 34

4 After-Sales Service

This material is under a two year warranty and should be returned to our stores in the event of any defects.

For any repairs, adjustments or spare parts, please contact:

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SA capital 3 233 762 € - Siren R.C.S. B 387 901 044 - Siret 387 901 04400017

Direct connection for technical support

A team of experts at your disposal from Monday to Friday (opening hours)

- You're looking for technical information ?
- You wish advice for use ?
- You need an urgent diagnosis ?

We take in charge your request immediately to provide you with the right answers regarding your activity field : Biology, Physics, Chemistry, Technology .

Free service *

+ 33 (0)2 32 29 40 50**

* Call cost not included

** Only for call from foreign countries